

PRIMO[®]



The perfect hydrocolloid tapered edge for
outstanding adhesion and security

OakMed 

Giving you more choice.

The unique features and benefits of the new Primo® range are integral to both closed and drainable pouches.



High tack adhesive

Providing state-of-the-art, infinitely superior adhesion.




Split fabric front cover

On **all** pouches, for accurate positioning and easy observation.



New, thicker, superior soft fabric cover

Won't cling to clothing and will towel dry after showering or bathing.



New anatomically designed shape

To fit body contours more closely.

PRIMO® 


The perfect hydrocolloid tapered edge adhesive

Hydrocolloid Adhesive

Hypo-allergenic, gentle on skin and non-irritant.



HYDROCOLLOID®




Unique tapered edge border

Will mould, shape and sculpt over body contours, reducing the need for flange extension products.



Two types of filter systems

Designed individually to meet the needs of colostomies and/or ileostomies.



Generous adhesive area

Gives security and confidence, even with larger stomas. When pouch is cut to the maximum 60mm, there is still 3cm of adhesive area left.



SafeSeal hygienic closure

Exclusive to our drainable pouches. Washable and can be tucked inside soft cover.



The perfect hydrocolloid tapered edge adhesive

Tapered edge hydrocolloid adhesive

The Primo[®] hydrocolloid adhesive provides: state of the art, infinitely superior, high tack adhesion, which in turn gives a faultless secure fit over body contours. The quality of the adhesion is unequalled, perfect in its composition.

The Primo[®] hydrocolloid adhesive will mould and stick instantly to the skin but can be removed and re-positioned whilst still maintaining the high-tack adhesion.

With its unique tapered edge integral border, the Primo[®] hydrocolloid adhesive undoubtedly reduces the need for flange extension products which makes pouch changes simpler and more cost effective.

Because of its unique mould, shape and sculpt features, the Primo[®] hydrocolloid will allow for extended wear time reducing the number of pouch changes required.

The Primo[®] hydrocolloid is hypo-allergenic, gentle on the skin, non-irritant and is easily removed leaving little or no residue behind.



Superior adhesion and tack



HYDROCOLLOID[®]

PRIMO®

The perfect hydrocolloid tapered edge adhesive

Brand new anatomically shaped closed and drainable pouches

Designed with comfort and confidence in mind

Primo® pouches have been designed with a unique anatomical shape to fit around body contours faultlessly, providing the wearer with increased protection and confidence.

- Improves comfort and fit when wearing the pouch
- Allows the pouch to tuck neatly into undergarments
- Creates an altogether more discreet experience



Large adhesive area. (3cm) still left, even when cut to 60mm.

Anatomically shaped




The perfect hydrocolloid tapered edge adhesive

Filters

Primo® drainable and closed pouches feature carbon filters which eliminate odour, helping to reduce ballooning and pancaking.


Primo® drainable pouch

Features the new Moisture Protected Composition (MPC) filter. This filter, along with its covering Gortex™ film, is specifically designed with an absorbent membrane to be liquid tight. This enables the filter to cope with the increased fluid output from an ileostomy.



Primo® closed pouch

A filter with a covering Gortex™ film, similar in composition to the MPC, but with less layers, which allows gas to pass through but not liquid, reducing spotting onto clothing.





The perfect hydrocolloid tapered edge adhesive

Split fabric cover



The Primo® range now features the NEW thicker, softer, more luxurious and comfortable SPLIT fabric cover.

This fabric cover is quiet, won't cling to clothing and will towel dry after showering or bathing.

The fabric is split over the stoma area to assist with application of the pouch and also allows viewing of the stoma area if required.

Luxurious fabric cover

- Won't cling to clothing
- Will towel dry after showering or bathing
- Is quiet, comfortable and discreet



Luxurious fabric cover

PRIMO®

Hydrocolloid Primo® Pouches

Standard Closed Pouches	Packs of 30
Hole Size	Opaque / Window
Cut to fit 20mm to 60mm	HP5C20
25mm	HP5C25
30mm	HP5C30
35mm	HP5C35
40mm	HP5C40
45mm	HP5C45
50mm	HP5C50
55mm	HP5C55
60mm	HP5C60

Opaque, window, filter; soft cover both sides

Standard Drainable Pouches	Packs of 30
Hole Size	Opaque / Window
Cut to fit 20mm to 60mm	HP8D20
25mm	HP8D25
30mm	HP8D30
35mm	HP8D35
40mm	HP8D40
45mm	HP8D45
50mm	HP8D50
55mm	HP8D55
60mm	HP8D60

Opaque, window, filter; soft cover both sides, velcro

OakMed 
Giving you more choice.



Try our fast and secure delivery service.
FREE complimentary items.
For example toiletry bags, mattress covers and radar keys.


OakMedGold
Your home delivery service

Oakmed Limited, Unit B, The I O Centre, Barn Way,
Lodge Farm Industrial Estate, Northampton, NN5 7UW

Freephone: 0800 592 786 | Main: 01604 586 529 | Fax: 01604 584 784 www.oakmed.co.uk

