



# PRODUCT GUIDE 2025

Edition 2




[www.oakmed.co.uk](http://www.oakmed.co.uk)

# Oakmed Labels Explained

Please use this chart to help you identify individual products within the Oakmed range

- All Oakmed products are colour coded and have a colour coded logo.
- You will see these colours and logos throughout our labelling, product guides, marketing and website.
- Some of the labels may have two colours if they come into two categories, eg. Soft Convex®/Alginate.
- All stock sent into hospitals are half stock boxes and state that on the label.

NB. Your Oakmed Representative is happy to check Oakmed stock for you when they visit.



## Reliability has a new colour

Discover our new range of black bags. Now available in our Soft Convex and Low Profile Soft Convex Drainable ranges; in both our gentle hydrocolloid tapered edge and naturally soothing alginate options.

Our new black bags will also be available in our equivalent closed ranges later in 2025.



**Ultimate softness**

Soft, smooth textile cover



**Hydrophobic surface**

Water-repellent



**Blackout weave**

100% Opaque



**Designed for comfort**

Asymmetrical shape

Find our new black bags on pages 14-16, 18, 20-22 and 26-28

## EXTENDING OUR RANGE OF DEPENDABLE SILICONE ACCESSORIES

Trio Sil2™ silicone accessories are now available from Oakmed and Goldcare

This range of accessories uses unique Sil2™ Breathable Silicone Technology, by Trio, laying the foundations for healthy stoma skin that is continually protected.

The range includes seals, sprays, wipes, gels, and flange extenders, designed to support essential skincare needs for ostomates, no matter how frequently they need to change their bag.

**NEW  
for  
2025**



No change, just confidence.

Find our new range of accessories on pages 49-52

## Starter Pack

Clinic/Community

### Clinic/Community closed bag starter packs

in all standard and maxi cut to fit closed bags

These packs are perfect for helping colostomists find the bag that suits them best.

Each pack contains 10 bags of the selected product — giving colostomists plenty of time to assess whether the bag meets their needs.

Full instructions for use are included to provide guidance on application and removal, along with a sizing guide to ensure the bag wafer is cut to the correct size.

The cover card details all the starter packs available and may be used as a reference tool to choose an alternative bag, if the initial one chosen does not suit the patient.




Ask your Oakmed Territory Manager for more details


The Oakmed bags are developed and manufactured in Berlin using modern technologies and high-quality materials. They are characterised by innovative and unique properties that meet the medical needs of healthcare professionals and patients.

## Bag shapes

Symmetric



Asymmetric




## Bag sizes

Drainable

Closed



## Closures



## Cover Material



Non-woven cover



Textile cover



NEW



## "Window"-Technology



All of the Oakmed bags feature a window for easy positioning and observation. Our clear bags without cover/window feature a crystal clear film, again for easy observation.

Filter




Activated carbon filter – Colostomy

MPC filter – Ileostomy

MPC filter technology (Moisture Protect Composition) = additional protective layer that prevents liquid absorption.

## Convexities



Low Profile Soft Convex

Soft Convex

# Contents



## Primo®

The perfect Hydrocolloid tapered edge for outstanding adhesion and security.

Page 10



## Soft Convex®

7mm convexity.

Available in tapered edge Hydrocolloid and Alginate adhesives. Also high output.

Page 13



## Low Profile Soft Convex®

3.5mm convexity.

Available in tapered edge Hydrocolloid and Alginate adhesives.

Page 17



## Alginate

An adhesive which contains Alginate, known for its moisture absorbing and soothing properties.

Page 23



## Urostomy

Available in tap and soft bung with Hydrocolloid, Microskin, Soft Convex and Low Profile Soft Convex adhesives.

Page 29



## Microskin®

The lightest, thinnest and most flexible wafer available.

Page 39



## Hydrocolloid

A soft, flexible, hypoallergenic adhesive.

Page 41



## Wound Care

Wound and fistula drainage appliances in a range of sizes.

Page 44



## Accessories

A comprehensive range of stoma and related accessories.

Page 48

Primo

Soft Convex

Low Profile Soft Convex

Alginate

Urostomy

Microskin

Hydrocolloid

Wound Care

Accessories

# PRIMO®

Designed by patients and Healthcare Professionals, the unique features and benefits of our Primo® bags offer you the perfect fit with outstanding adhesion and assured security.

These features are integral to both our drainable and closed bags.

- Oakmed's most popular range, Primo® has the perfect Hydrocolloid tapered edge baseplate for outstanding adhesion and security.
- Anatomically shaped.
- Infinitely superior, high tack adhesion.
- Leaves little or no residue when removed.
- Unique tapered edge, integral border.
- Gentle Hydrocolloid.
- Large adhesive area even when cut to maximum.
- Carbon filters designed individually for closed and drainable bag requirements.
- Improved soft touch closure on drainable bags.
- Soft non-cling fabric cover with viewing window.

**Our most popular range with mini, midi, standard and maxi sizes available. Soft touch, easy to use closure on drainable bags.**



The perfect fit with outstanding adhesion.



## Closed Mini Hydrocolloid Primo

HP11C10 Window, opaque, filter, soft cover both sides.

Pack Size: 30

Hole Size	Window
Cut to fit 10 to 50mm	HP11C10



## Closed Standard Hydrocolloid Primo

HP5C20 - HP5C60 Window, opaque, filter, soft cover both sides.

HP5C20C Clear, filter, soft cover one side.

Pack Size: 30

Hole Size	Window	Clear
Cut to fit 20 to 60mm	HP5C20	HP5C20C
Pre-cut to 25mm	HP5C25	
Pre-cut to 30mm	HP5C30	
Pre-cut to 35mm	HP5C35	
Pre-cut to 40mm	HP5C40	
Pre-cut to 45mm	HP5C45	
Pre-cut to 50mm	HP5C50	
Pre-cut to 55mm	HP5C55	
Pre-cut to 60mm	HP5C60	



## Closed Maxi Hydrocolloid Primo

HP8C20 Window, opaque, filter, soft cover both sides.

HP7C20C Clear, filter, soft cover one side.

Pack Size: 30

Hole Size	Window
Cut to fit 20 to 90mm	HP8C20

The perfect fit with outstanding adhesion.



### Drainable Midi Hydrocolloid Primo

HP8D20 - HP8D60 Window, opaque, filter, soft cover both sides, soft touch closure.

HP8D20C Clear, filter, soft cover one side, soft touch closure.

Pack Size: 30

Hole Size	Window	Clear
Cut to fit 20 to 60mm	HP8D20	HP8D20C
Pre-cut to 25mm	HP8D25	
Pre-cut to 30mm	HP8D30	
Pre-cut to 35mm	HP8D35	
Pre-cut to 40mm	HP8D40	
Pre-cut to 45mm	HP8D45	
Pre-cut to 50mm	HP8D50	
Pre-cut to 55mm	HP8D55	
Pre-cut to 60mm	HP8D60	



### Drainable Standard Hydrocolloid Primo

HP10D20 - HP10D60 Window, opaque, filter, soft cover both sides, soft touch closure.

HP10D20C Clear, filter, soft cover one side, soft touch closure.

Pack Size: 30


Hole Size	Window	Clear
Cut to fit 10 to 50mm	HP10D20	HP10D20C
Pre-cut to 25mm	HP10D25	
Pre-cut to 30mm	HP10D30	
Pre-cut to 35mm	HP10D35	
Pre-cut to 40mm	HP10D40	
Pre-cut to 45mm	HP10D45	
Pre-cut to 50mm	HP10D50	
Pre-cut to 55mm	HP10D55	
Pre-cut to 60mm	HP10D60	

## 7mm convexity

Naturally the flexible choice for maximum security and added confidence, our Soft Convex® range of bags fit the contours of your body for more comfortable wear. In both drainable and closed, tapered edge Hydrocolloid and Alginate adhesives, our range also offers a high output bag.

Unique features:

- 7mm depth of convexity.
- Maintains its convexity whilst in use.
- A comprehensive range of cut to fit and pre-cut sizes.
- High tack adhesion for advanced performance.
- Skin friendly design, hypoallergenic adhesives.
- Efficient odour eliminating filter system.




**Our Drainable Soft Convex bags are now available in black.**  
**Ultimate softness with a black out woven fabric that is water repellent.**  
**In both alginate and hydrocolloid for trusted skin protection.**





Naturally the flexible choice.



COMING SOON

### Closed Standard Hydrocolloid Soft Convex

SC25-0410KBK - SC57-0420K Window, opaque filter, soft cover both sides.

Pack Size: 20

Hole Size	Plateau Size	Window
<b>Cut to fit 10 to 24mm</b>	<b>25mm</b>	<b>SC25-0410KBK</b>
Cut to fit 20 to 24mm	25mm	SC25-0420K
<b>Cut to fit 15 to 34mm</b>	<b>35mm</b>	<b>SC35-0415KBK</b>
Cut to fit 20 to 34mm	35mm	SC35-0420K
Pre-cut to 25mm	35mm	SC35-0425K
Pre-cut to 30mm	35mm	SC35-0430K
<b>Cut to fit 20 to 44mm</b>	<b>45mm</b>	<b>SC45-0420KBK</b>
Cut to fit 20 to 44mm	45mm	SC45-0420K
Pre-cut to 35mm	45mm	SC45-0435K
Pre-cut to 40mm	45mm	SC45-0440K

Pack Size: 10

Cut to fit 20 to 57mm	57mm	SC57-0420K
-----------------------	------	------------

### Drainable Midi Hydrocolloid Soft Convex

SC25-4315KV - SC45-4340KV Window, opaque, filter, soft cover both sides, soft touch closure.

SC25-4315KVC - SC45-4315KVC Clear, filter, soft cover one side, soft touch closure.


Pack Size: 20

Hole Size	Plateau Size	Window
Cut to fit 15 to 24mm	25mm	SC25-4315KV
Pre-cut to 20mm	25mm	SC25-4320KV
Cut to fit 15 to 34mm	35mm	SC35-4315KV
Pre-cut to 25mm	35mm	SC35-4325KV
Pre-cut to 30mm	35mm	SC35-4330KV
Cut to fit 15 to 44mm	45mm	SC45-4315KV
Pre-cut to 35mm	45mm	SC45-4335KV
Pre-cut to 40mm	45mm	SC45-4340KV

Hole Size	Plateau Size	Clear
Cut to fit 15 to 24mm	25mm	SC25-4315KVC
Cut to fit 15 to 34mm	35mm	SC35-4315KVC
Cut to fit 15 to 44mm	45mm	SC45-4315KVC



Naturally the flexible choice.



### Drainable Standard Hydrocolloid Soft Convex

SC25-4110KBK - SC57-4120KV Window, opaque, filter, soft cover both sides, soft touch closure.

Pack Size: 20

Hole Size	Plateau Size	Window
<b>Cut to fit 10 to 24mm</b>	<b>25mm</b>	<b>SC25-4110KBK</b>
Cut to fit 20 to 24mm	25mm	SC25-4120KV
<b>Cut to fit 15 to 34mm</b>	<b>35mm</b>	<b>SC35-4115KBK</b>
Cut to fit 20 to 34mm	35mm	SC35-4120KV
Pre-cut to 25mm	35mm	SC35-4125KV
Pre-cut to 30mm	35mm	SC35-4130KV
<b>Cut to fit 15 to 44mm</b>	<b>45mm</b>	<b>SC45-4115KBK</b>
Cut to fit 20 to 44mm	45mm	SC45-4120KV
Pre-cut to 35mm	45mm	SC45-4135KV
Pre-cut to 40mm	45mm	SC45-4140KV

Pack Size: 10

Cut to fit 20 to 57mm	57mm	SC57-4120KV
-----------------------	------	-------------

### Drainable Standard Hydrocolloid High Output Soft Convex

SCHV25-4115K - SCHV45-4115K Window, opaque, filter, soft cover both sides, wide bung closure.

Pack Size: 10


Hole Size	Plateau Size	Window
Cut to fit 15 to 24mm	25mm	SCHV25-4115K
Cut to fit 15 to 34mm	35mm	SCHV35-4115K
Cut to fit 15 to 44mm	45mm	SCHV45-4115K



All bag types in **bold** are in black cover

All bag types in **bold** are in black cover

Naturally the flexible choice.



COMING SOON

### Closed Standard Alginate Soft Convex

SCA25-0410KBK - SCA45-0440K Window, opaque filter, soft cover both sides.

Pack Size: 20

Hole Size	Plateau Size	Window
<b>Cut to fit 10 to 24mm</b>	<b>25mm</b>	<b>SCA25-0410KBK</b>
Cut to fit 20 to 24mm	25mm	SCA25-0420K
<b>Cut to fit 15 to 34mm</b>	<b>35mm</b>	<b>SCA35-0415KBK</b>
Cut to fit 20 to 34mm	35mm	SCA35-0420K
Pre-cut to 25mm	35mm	SCA35-0425K
Pre-cut to 30mm	35mm	SCA35-0430K
<b>Cut to fit 20 to 44mm</b>	<b>45mm</b>	<b>SCA45-0420KBK</b>
Cut to fit 20 to 44mm	45mm	SCA45-0420K
Pre-cut to 35mm	45mm	SCA45-0435K
Pre-cut to 40mm	45mm	SCA45-0440K

### Drainable Standard Alginate Soft Convex

SCA25-4110KVBK - SCA45-4140KV Window, opaque, filter, soft cover both sides, soft touch closure.

Pack Size: 20

Hole Size	Plateau Size	Window
<b>Cut to fit 10 to 24mm</b>	<b>25mm</b>	<b>SCA25-4110KVBK</b>
Cut to fit 20 to 24mm	25mm	SCA25-4120KV
<b>Cut to fit 15 to 34mm</b>	<b>35mm</b>	<b>SCA35-4115KVBK</b>
Cut to fit 20 to 34mm	35mm	SCA35-4120KV
Pre-cut to 25mm	35mm	SCA35-4125KV
Pre-cut to 30mm	35mm	SCA35-4130KV
<b>Cut to fit 15 to 44mm</b>	<b>45mm</b>	<b>SCA45-4115KVBK</b>
Cut to fit 20 to 44mm	45mm	SCA45-4120KV
Pre-cut to 35mm	45mm	SCA45-4135KV
Pre-cut to 40mm	45mm	SCA45-4140KV

All bag types in **bold** are in black cover

low profile  
soft convex®

Low Profile Soft Convex


### 3.5mm convexity

Offering the same superior flexibility and comfortable wear, our Low Profile Soft Convex® range gives you 3.5mm of convexity.

Available in both drainable and closed, tapered edge Hydrocolloid and Alginate adhesives.

Unique features:

- 3.5mm depth of convexity.
- Extremely cost effective product.
- Soft and flexible design, works like a flat bag with a seal.
- A comprehensive range of cut to fit and pre-cut sizes.
- High tack adhesion for advanced performance.
- Skin friendly design, hypoallergenic adhesives.
- Efficient odour eliminating filter system.



**NEW**  
Black bags


Our Drainable Low Profile Soft Convex bags are now available in black. Ultimate softness with a black out woven fabric that is water repellent. In both alginate and hydrocolloid for trusted skin protection.



CLOSED  
BLACK BAGS  
COMING  
SOON



Naturally the flexible choice.



COMING  
SOON



### Closed Standard Hydrocolloid Low Profile Soft Convex


LPSC28-0410KBK - LPSC47-0445K Window, opaque, filter, soft cover both sides.

Pack Size: 20

Hole Size	Plateau Size	Window
<b>Cut to fit 10 to 28mm</b>	<b>28mm</b>	<b>LPSC28-0410KBK</b>
Cut to fit 20 to 28mm	28mm	LPSC28-0420K
Pre-cut to 25mm	28mm	LPSC28-0425K
<b>Cut to fit 15 to 40mm</b>	<b>40mm</b>	<b>LPSC40-0415KBK</b>
Cut to fit 20 to 40mm	40mm	LPSC40-0420K
Pre-cut to 30mm	40mm	LPSC40-0430K
Pre-cut to 35mm	40mm	LPSC40-0435K
<b>Cut to fit 20 to 47mm</b>	<b>47mm</b>	<b>LPSC47-0420KBK</b>
Cut to fit 20 to 47mm	47mm	LPSC47-0420K
Pre-cut to 40mm	47mm	LPSC47-0440K
Pre-cut to 45mm	47mm	LPSC47-0445K



Naturally the flexible choice.



### Drainable Midi Hydrocolloid Low Profile Soft Convex

LPSC28-4315KV - LPSC47-4345KV Window, opaque, filter, soft cover both sides, soft touch closure.

LPSC28-4315KVC - LPSC47-4315KVC Clear, filter, soft cover one side, soft touch closure.


Pack Size: 20

Hole Size	Plateau Size	Window	Clear
Cut to fit 15 to 28mm	28mm	LPSC28-4315KV	LPSC28-4315KVC
Pre-cut to 20mm	28mm	LPSC28-4320KV	
Pre-cut to 25mm	28mm	LPSC28-4325KV	
Cut to fit 15 to 40mm	40mm	LPSC40-4315KV	LPSC40-4315KVC
Pre-cut to 30mm	40mm	LPSC40-4330KV	
Pre-cut to 35mm	40mm	LPSC40-4335KV	
Cut to fit 15 to 47mm	47mm	LPSC47-4315KV	LPSC47-4315KVC
Pre-cut to 40mm	47mm	LPSC47-4340KV	
Pre-cut to 45mm	47mm	LPSC47-4345KV	

All bag types in **bold** are in black cover



Naturally the flexible choice.



### Drainable Standard Hydrocolloid Low Profile Soft Convex

LPSC28-4110KVBK - LPSC47-4145KV Window, opaque, filter, soft cover both sides, soft touch closure.

LPSC28-4115KVC - LPSC47-4115KVC Clear, filter, soft cover one side, soft touch closure.



Pack Size: 20

Hole Size	Plateau Size	Window	Clear
<b>Cut to fit 10 to 28mm</b>	<b>28mm</b>	<b>LPSC28-4110KVBK</b>	
Cut to fit 15 to 28mm	28mm	LPSC28-4115KV	LPSC28-4115KVC
Pre-cut to 20mm	28mm	LPSC28-4120KV	
Pre-cut to 25mm	28mm	LPSC28-4125KV	
<b>Cut to fit 15 to 40mm</b>	<b>40mm</b>	<b>LPSC40-4115KVBK</b>	
Cut to fit 15 to 40mm	40mm	LPSC40-4115KV	LPSC40-4115KVC
Pre-cut to 30mm	40mm	LPSC40-4130KV	
Pre-cut to 35mm	40mm	LPSC40-4135KV	
<b>Cut to fit 15 to 47mm</b>	<b>47mm</b>	<b>LPSC47-4115KVBK</b>	
Cut to fit 15 to 47mm	47mm	LPSC47-4115KV	LPSC47-4115KVC
Pre-cut to 40mm	47mm	LPSC47-4140KV	
Pre-cut to 45mm	47mm	LPSC47-4145KV	

All bag types in **bold** are in black cover



Naturally the flexible choice.



### Closed Standard Alginate Low Profile Soft Convex

LPSCA28-0410KBK - LPSCA47-0445K Window, opaque, filter, soft cover both sides.


Pack Size: 20

Hole Size	Plateau Size	Window
<b>Cut to fit 10 to 28mm</b>	<b>28mm</b>	<b>LPSCA28-0410KBK</b>
Cut to fit 20 to 28mm	28mm	LPSCA28-0420K
Pre-cut to 25mm	28mm	LPSCA28-0425K
<b>Cut to fit 15 to 40mm</b>	<b>40mm</b>	<b>LPSCA40-0415KBK</b>
Cut to fit 20 to 40mm	40mm	LPSCA40-0420K
Pre-cut to 30mm	40mm	LPSCA40-0430K
Pre-cut to 35mm	40mm	LPSCA40-0435K
<b>Cut to fit 20 to 47mm</b>	<b>47mm</b>	<b>LPSCA47-0420KBK</b>
Cut to fit 20 to 47mm	47mm	LPSCA47-0420K
Pre-cut to 40mm	47mm	LPSCA47-0440K
Pre-cut to 45mm	47mm	LPSCA47-0445K

All bag types in **bold** are in black cover



Naturally the flexible choice.



### Drainable Standard Alginate Low Profile Soft Convex

LPSCA28-4110KVBK - LPSCA47-4145KV Window, opaque, filter, soft cover both sides, soft touch closure.

Pack Size: 20

Hole Size	Plateau Size	Window
<b>Cut to fit 10 to 28mm</b>	<b>28mm</b>	<b>LPSCA28-4110KVBK</b>
Cut to fit 15 to 28mm	28mm	LPSCA28-4115KV
Pre-cut to 20mm	28mm	LPSCA28-4120KV
Pre-cut to 25mm	28mm	LPSCA28-4125KV
<b>Cut to fit 15 to 40mm</b>	<b>40mm</b>	<b>LPSCA40-4115KVBK</b>
Cut to fit 15 to 40mm	40mm	LPSCA40-4115KV
Pre-cut to 30mm	40mm	LPSCA40-4130KV
Pre-cut to 35mm	40mm	LPSCA40-4135KV
<b>Cut to fit 15 to 47mm</b>	<b>47mm</b>	<b>LPSCA47-4115KVBK</b>
Cut to fit 15 to 47mm	47mm	LPSCA47-4115KV
Pre-cut to 40mm	47mm	LPSCA47-4140KV
Pre-cut to 45mm	47mm	LPSCA47-4145KV

All bag types in **bold** are in black cover



### A smart choice for problem skin

Alginate is a natural ingredient derived from seaweed, known for its moisture absorbing and soothing properties. Our Alginate baseplate is infused with 15% Alginate and may help to promote skin health around the stoma.

Many ostomates face skin issues from moisture build-up, leaks, friction or allergies to other products, leading to painful rashes and inflammation, affecting comfort and quality of life. Our alginate range is available in closed and drainable options.

#### Benefits of Oakmed's Alginate range

- Hydrophilic – may absorb and lock away moisture for a dry and comfortable wear.
- Alginate, a natural ingredient derived from seaweed.
- Anatomical shape for comfort and security.
- Warm up to wake up - place under your arm or on a radiator to wake up the soothing properties of the wafer.
- Secure and strong closure for reliability and security, giving you peace of mind.

Oakmed Alginate adhesives contain not just a trace element but **15% Alginate**.





With **15% Alginate**, a smart choice for problem skin



### Closed Mini Alginate

AHP11C10 Window, opaque, filter, soft cover both sides.

Pack Size: 30

Hole Size	Window
Cut to fit 10 to 50mm	AHP11C10



With **15% Alginate**, a smart choice for problem skin



### Drainable Midi Alginate

AV-4320K Window, filter, soft cover both sides, soft touch closure.

Pack Size: 30

Hole Size	Window Midi
Cut to fit 20 to 60mm	AV-4320K



### Closed Standard Alginate

AHP5C20 - AHP5C60 Window, opaque, filter, soft cover both sides.

Pack Size: 30

Hole Size	Window
Cut to fit 20 to 60mm	AHP5C20
Pre-cut to 25mm	AHP5C25
Pre-cut to 30mm	AHP5C30
Pre-cut to 35mm	AHP5C35
Pre-cut to 40mm	AHP5C40
Pre-cut to 45mm	AHP5C45
Pre-cut to 50mm	AHP5C50
Pre-cut to 55mm	AHP5C55
Pre-cut to 60mm	AHP5C60



### Drainable Standard Alginate

AV-4120K - AV-4160K Window, filter, soft cover both sides, soft touch closure.

Pack Size: 30

Hole Size	Window
Cut to fit 20 to 60mm	AV-4120K
Pre-cut to 25mm	AV-4125K
Pre-cut to 30mm	AV-4130K
Pre-cut to 35mm	AV-4135K
Pre-cut to 40mm	AV-4140K
Pre-cut to 45mm	AV-4145K
Pre-cut to 50mm	AV-4150K
Pre-cut to 55mm	AV-4155K
Pre-cut to 60mm	AV-4160K



### Closed Maxi Alginate

AHP8C20 Window, opaque, filter, soft cover both sides.


AHP7C20C Clear, soft cover one side.

Pack Size: 30

Hole Size	Window
Cut to fit 20 to 90mm	AHP8C20



Naturally the flexible choice.



COMING SOON

### Closed Standard Alginate Soft Convex

SCA25-0410KBK - SCA45-0440K Window, opaque filter, soft cover both sides.

Pack Size: 20

Hole Size	Plateau Size	Window
<b>Cut to fit 10 to 24mm</b>	<b>25mm</b>	<b>SCA25-0410KBK</b>
Cut to fit 20 to 24mm	25mm	SCA25-0420K
<b>Cut to fit 15 to 34mm</b>	<b>35mm</b>	<b>SCA35-0415KBK</b>
Cut to fit 20 to 34mm	35mm	SCA35-0420K
Pre-cut to 25mm	35mm	SCA35-0425K
Pre-cut to 30mm	35mm	SCA35-0430K
<b>Cut to fit 20 to 44mm</b>	<b>45mm</b>	<b>SCA45-0420KBK</b>
Cut to fit 20 to 44mm	45mm	SCA45-0420K
Pre-cut to 35mm	45mm	SCA45-0435K
Pre-cut to 40mm	45mm	SCA45-0440K

### Drainable Standard Alginate Soft Convex

SCA25-4110KBK - SCA45-4140KV Window, opaque, filter, soft cover both sides, soft touch closure.

Pack Size: 20




Hole Size	Plateau Size	Window
<b>Cut to fit 10 to 24mm</b>	<b>25mm</b>	<b>SCA25-4110KBK</b>
Cut to fit 20 to 24mm	25mm	SCA25-4120KV
<b>Cut to fit 15 to 34mm</b>	<b>35mm</b>	<b>SCA35-4115KBK</b>
Cut to fit 20 to 34mm	35mm	SCA35-4120KV
Pre-cut to 25mm	35mm	SCA35-4125KV
Pre-cut to 30mm	35mm	SCA35-4130KV
<b>Cut to fit 15 to 44mm</b>	<b>45mm</b>	<b>SCA45-4115KBK</b>
Cut to fit 20 to 44mm	45mm	SCA45-4120KV
Pre-cut to 35mm	45mm	SCA45-4135KV
Pre-cut to 40mm	45mm	SCA45-4140KV

All bag types in **bold** are in black cover



Naturally the flexible choice.



COMING SOON



### Closed Standard Alginate Low Profile Soft Convex

LPSCA28-0410KBK - LPSCA47-0445K Window, opaque, filter, soft cover both sides.


Pack Size: 20

Hole Size	Plateau Size	Window
<b>Cut to fit 10 to 28mm</b>	<b>28mm</b>	<b>LPSCA28-0410KBK</b>
Cut to fit 20 to 28mm	28mm	LPSCA28-0420K
Pre-cut to 25mm	28mm	LPSCA28-0425K
<b>Cut to fit 15 to 40mm</b>	<b>40mm</b>	<b>LPSCA40-0415KBK</b>
Cut to fit 20 to 40mm	40mm	LPSCA40-0420K
Pre-cut to 30mm	40mm	LPSCA40-0430K
Pre-cut to 35mm	40mm	LPSCA40-0435K
<b>Cut to fit 20 to 47mm</b>	<b>47mm</b>	<b>LPSCA47-0420KBK</b>
Cut to fit 20 to 47mm	47mm	LPSCA47-0420K
Pre-cut to 40mm	47mm	LPSCA47-0440K
Pre-cut to 45mm	47mm	LPSCA47-0445K

All bag types in **bold** are in black cover



Naturally the flexible choice.



### Drainable Standard Alginate Low Profile Soft Convex

LPSCA28-4110KVBK - LPSCA47-4145KV Window, opaque, filter, soft cover both sides, soft touch closure.

Pack Size: 20

Hole Size	Plateau Size	Window
<b>Cut to fit 10 to 28mm</b>	<b>28mm</b>	<b>LPSCA28-4110KVBK</b>
Cut to fit 15 to 28mm	28mm	LPSCA28-4115KV
Pre-cut to 20mm	28mm	LPSCA28-4120KV
Pre-cut to 25mm	28mm	LPSCA28-4125KV
<b>Cut to fit 15 to 40mm</b>	<b>40mm</b>	<b>LPSCA40-4115KVBK</b>
Cut to fit 15 to 40mm	40mm	LPSCA40-4115KV
Pre-cut to 30mm	40mm	LPSCA40-4130KV
Pre-cut to 35mm	40mm	LPSCA40-4135KV
<b>Cut to fit 15 to 47mm</b>	<b>47mm</b>	<b>LPSCA47-4115KVBK</b>
Cut to fit 15 to 47mm	47mm	LPSCA47-4115KV
Pre-cut to 40mm	47mm	LPSCA47-4140KV
Pre-cut to 45mm	47mm	LPSCA47-4145KV

All bag types in **bold** are in black cover



Our extended Urostomy range is shaped for maximum capacity with minimum bulk. The unique anti-twist shape allows for easier emptying. Now available in soft bung or new tap closure for added confidence and convenience.

- Available in a choice of adhesives.
- Soft fabric cover, quiet, non-cling, comfortable.
- Split fabric cover for easy positioning.
- Non-return valve, stops urine back flowing.
- Soft bung closure, safe from leaks.
- Anti-twist shape.
- Drip cap on tap closure.

Ask about our Clinic/Community packs, the easy way for patients to find the right bag for them.






Maximum capacity, minimum bulk for increased confidence and comfort.



Maximum capacity, minimum bulk for increased confidence and comfort.



### Urostomy Standard Hydrocolloid

GC1500 Window, opaque, soft cover both sides, soft bung closure, capacity 420ml.

Pack Size: 20

Hole Size	Window
✦ Cut to fit 15 to 60mm	GC1500



### Urostomy Standard Hydrocolloid

GC1500T Window, opaque, soft cover both sides, tap closure, capacity 420ml.

Pack Size: 20

Hole Size	Window
✦ Cut to fit 15 to 60mm	GC1500T



### Urostomy Standard Hydrocolloid

GC1500C Clear, soft cover one side, soft bung closure, capacity 420ml.

Pack Size: 20

Hole Size	Clear
✦ Cut to fit 15 to 60mm	GC1500C



### Urostomy Standard Hydrocolloid

GC1500CT Clear, soft cover both sides, tap closure, capacity 420ml.

Pack Size: 20

Hole Size	Clear
✦ Cut to fit 15 to 60mm	GC1500CT

✦ Available as a Clinic/Community pack. Call us to order.

✦ Available as a Clinic/Community pack. Call us to order.



Maximum capacity, minimum bulk, with 7mm convexity, for increased confidence and comfort.



### Urostomy Standard Hydrocolloid Soft Convex

SCU25-0015 - SCU45-0040 Window, opaque, soft cover both sides, soft bung closure, capacity 420ml.

Pack Size: 20



Hole Size	Plateau Size	Window
* Cut to fit 15 to 24mm	25mm	SCU25-0015
Pre-cut to 20mm	25mm	SCU25-0020
* Cut to fit 15 to 34mm	35mm	SCU35-0015
Pre-cut to 25mm	35mm	SCU35-0025
Pre-cut to 30mm	35mm	SCU35-0030
* Cut to fit 15 to 44mm	45mm	SCU45-0015
Pre-cut to 35mm	45mm	SCU45-0035
Pre-cut to 40mm	45mm	SCU45-0040

### Urostomy Standard Hydrocolloid Soft Convex

SCU25-0015C - SCU45-0015C Clear, soft cover one side, soft bung closure, capacity 420ml.

Pack Size: 20




Hole Size	Plateau Size	Clear
* Cut to fit 15 to 24mm	25mm	SCU25-0015C
* Cut to fit 15 to 34mm	35mm	SCU35-0015C
* Cut to fit 15 to 44mm	45mm	SCU45-0015C

\* Available as a Clinic/Community pack. Call us to order.



Maximum capacity, minimum bulk, with 7mm convexity, for increased confidence and comfort.



### Urostomy Standard Hydrocolloid Soft Convex

SCU25-0015T - SCU45-0040T Window, opaque, soft cover both sides, tap closure, capacity 420ml.

Pack Size: 20



Hole Size	Plateau Size	Window
* Cut to fit 15 to 24mm	25mm	SCU25-0015T
Pre-cut to 20mm	25mm	SCU25-0020T
* Cut to fit 15 to 34mm	35mm	SCU35-0015T
Pre-cut to 25mm	35mm	SCU35-0025T
Pre-cut to 30mm	35mm	SCU35-0030T
* Cut to fit 15 to 44mm	45mm	SCU45-0015T
Pre-cut to 35mm	45mm	SCU45-0035T
Pre-cut to 40mm	45mm	SCU45-0040T

### Urostomy Standard Hydrocolloid Soft Convex

SCU25-0015CT - SCU45-0015CT Clear, soft cover one side, tap closure, capacity 420ml.

Pack Size: 20



Hole Size	Plateau Size	Clear
* Cut to fit 15 to 24mm	25mm	SCU25-0015CT
* Cut to fit 15 to 34mm	35mm	SCU35-0015CT
* Cut to fit 15 to 44mm	45mm	SCU45-0015CT

\* Available as a Clinic/Community pack. Call us to order.



Maximum capacity, minimum bulk, with 3.5mm convexity, for increased confidence and comfort.



**BUNG**

### Urostomy Standard Hydrocolloid Low Profile Soft Convex

LPSCU28-0015 - LPSCU47-0045 Window, opaque, soft cover both sides, bung closure, capacity 420ml.

Pack Size: 20



Hole Size	Plateau Size	Window
* Cut to fit 15 to 28mm	28mm	LPSCU28-0015
Pre-cut to 20mm	28mm	LPSCU28-0020
Pre-cut to 25mm	28mm	LPSCU28-0025
* Cut to fit 15 to 40mm	40mm	LPSCU40-0015
Pre-cut to 30mm	40mm	LPSCU40-0030
Pre-cut to 35mm	40mm	LPSCU40-0035
* Cut to fit 15 to 47mm	47mm	LPSCU47-0015
Pre-cut to 40mm	47mm	LPSCU47-0040
Pre-cut to 45mm	47mm	LPSCU47-0045

### Urostomy Standard Hydrocolloid Low Profile Soft Convex

LPSCU28-0015C - LPSCU47-0015C Clear, soft cover one side, soft bung closure, capacity 420ml.

Pack Size: 20



Hole Size	Plateau Size	Clear
* Cut to fit 15 to 28mm	28mm	LPSCU28-0015C
* Cut to fit 15 to 40mm	40mm	LPSCU40-0015C
* Cut to fit 15 to 47mm	47mm	LPSCU47-0015C


\* Available as a Clinic/Community pack. Call us to order.



Maximum capacity, minimum bulk, with 3.5mm convexity, for increased confidence and comfort.



**DRIP CAP**



**TAP**

### Urostomy Standard Hydrocolloid Low Profile Soft Convex

LPSCU28-0015T - LPSCU47-0045T Window, opaque, soft cover both sides, tap closure, capacity 420ml.

Pack Size: 20



Hole Size	Plateau Size	Window
* Cut to fit 15 to 28mm	28mm	LPSCU28-0015T
Pre-cut to 20mm	28mm	LPSCU28-0020T
Pre-cut to 25mm	28mm	LPSCU28-0025T
* Cut to fit 15 to 40mm	40mm	LPSCU40-0015T
Pre-cut to 30mm	40mm	LPSCU40-0030T
Pre-cut to 35mm	40mm	LPSCU40-0035T
* Cut to fit 15 to 47mm	47mm	LPSCU47-0015T
Pre-cut to 40mm	47mm	LPSCU47-0040T
Pre-cut to 45mm	47mm	LPSCU47-0045T

### Urostomy Standard Hydrocolloid Low Profile Soft Convex

LPSCU28-0015CT - LPSCU47-0015CT Clear, soft cover one side, tap closure, capacity 420ml.

Pack Size: 20



Hole Size	Plateau Size	Clear
* Cut to fit 15 to 28mm	28mm	LPSCU28-0015CT
* Cut to fit 15 to 40mm	40mm	LPSCU40-0015CT
* Cut to fit 15 to 47mm	47mm	LPSCU47-0015CT

\* Available as a Clinic/Community pack. Call us to order.

# Starter Pack

Clinic/community

## Clinic/Community starter packs in all cut to fit Urostomy products

We have a full range of clinic starter packs available in both new 'Twist and Safe' tap and soft bung closure options.

Ideal to help Urostomists choose the right bag for their needs. Healthcare Professionals can rest assured that the pack contains everything the patient needs to trial our bag options. Full instructions for use are included to provide guidance on how to use the connector with the night drainage bag.

Please see over the page for the options we offer.

Pack contains:

- 5 x Bags of chosen product
- 1 x Night drainage bag
- 2 x Connectors
- 1 x Belt
- 1 x IFU
- 1 x Product code information sheet



Ask your Oakmed Territory Manager for more details.

Please contact our customer service department on **0800 592 786** if you have any questions or need any help at all regarding this or any other Oakmed products.



### Tap

Hole Size	Window	Clear
✱ Cut to fit 15 to 60mm	GC1500T	-
✱ Cut to fit 15 to 60mm	-	GC1500CT
✱ Cut to fit 15 to 24mm	SCU25-0015T	-
✱ Cut to fit 15 to 34mm	SCU35-0015T	-
✱ Cut to fit 15 to 44mm	SCU45-0015T	-
✱ Cut to fit 15 to 24mm	-	SCU25-0015CT
✱ Cut to fit 15 to 34mm	-	SCU35-0015CT
✱ Cut to fit 15 to 44mm	-	SCU45-0015CT
✱ Cut to fit 15 to 28mm	LPSCU28-0015T	-
✱ Cut to fit 15 to 40mm	LPSCU40-0015T	-
✱ Cut to fit 15 to 47mm	LPSCU47-0015T	-
✱ Cut to fit 15 to 28mm	-	LPSCU28-0015CT
✱ Cut to fit 15 to 40mm	-	LPSCU40-0015CT
✱ Cut to fit 15 to 47mm	-	LPSCU47-0015CT

### Bung

Hole Size	Window	Clear
✱ Cut to fit 15 to 60mm	GC1500	-
✱ Cut to fit 15 to 60mm	-	GC1500C
✱ Cut to fit 15 to 24mm	SCU25-0015	-
✱ Cut to fit 15 to 34mm	SCU35-0015	-
✱ Cut to fit 15 to 44mm	SCU45-0015	-
✱ Cut to fit 15 to 24mm	-	SCU25-0015C
✱ Cut to fit 15 to 34mm	-	SCU35-0015C
✱ Cut to fit 15 to 44mm	-	SCU45-0015C
✱ Cut to fit 15 to 28mm	LPSCU28-0015	-
✱ Cut to fit 15 to 40mm	LPSCU40-0015	-
✱ Cut to fit 15 to 47mm	LPSCU47-0015	-
✱ Cut to fit 15 to 28mm	-	LPSCU28-0015C
✱ Cut to fit 15 to 40mm	-	LPSCU40-0015C
✱ Cut to fit 15 to 47mm	-	LPSCU47-0015C



Maximum capacity, minimum bulk for increased confidence and comfort.



### Urostomy Standard Microskin

MJ1300T - MJ3800T Window, opaque, soft cover both sides, tap closure.

Pack Size: 20

Hole Size	Window
Cut to fit 10 to 50mm	MJ1300T
Pre-cut to 16mm	MJ1600T
Pre-cut to 19mm	MJ1900T
Pre-cut to 22mm	MJ2200T
Pre-cut to 25mm	MJ2500T
Pre-cut to 32mm	MJ3200T
Pre-cut to 38mm	MJ3800T



### Urostomy Standard Microskin

MJ1300 - MJ3800 Window, opaque, soft cover both sides, soft bung closure.

Pack Size: 20

Hole Size	Window
Cut to fit 10 to 50mm	MJ1300
Pre-cut to 16mm	MJ1600
Pre-cut to 19mm	MJ1900
Pre-cut to 22mm	MJ2200
Pre-cut to 25mm	MJ2500
Pre-cut to 32mm	MJ3200
Pre-cut to 38mm	MJ3800

# microskin®

Unique to Oakmed, the Microskin® adhesive is thin, transparent and flexible, moulding easily to all body contours and creases, significantly reducing the risk of leakage.

Unique features:

- Transparent, thin and flexible design.
- Moulds easily and comfortably to body contours.
- Barely visible on skin when in place.
- Breathable and moisture vapour permeable for healthier skin.
- Durable for extended wear time.

Transparent and flexible, Microskin® is the ideal wafer for active people.



Transparent and flexible, reduces the risk of leakage.



### Closed Standard Microskin

M-0320K - M-0340K Window, filter, soft cover both sides.

Pack Size: 30

Hole Size	Window
Cut to fit 20 to 50mm	M-0320K
Pre-cut to 25mm	M-0325K
Pre-cut to 30mm	M-0330K
Pre-cut to 35mm	M-0335K
Pre-cut to 40mm	M-0340K



### Closed Plus Microskin

M-0420K - M-0440K Window, filter, soft cover both sides.

Pack Size: 30

Hole Size	Window
Cut to fit 20 to 50mm	M-0420K
Pre-cut to 25mm	M-0425K
Pre-cut to 30mm	M-0430K
Pre-cut to 35mm	M-0435K
Pre-cut to 40mm	M-0440K



### Drainable Standard Microskin

MV-4120K Window, filter, soft cover both sides, soft touch closure.

Pack Size: 30

Hole Size	Standard
Cut to fit 20 to 50mm	MV-4120K



## HYDROCOLLOID

Oakmed's Hydrocolloid has a long term reputation for minimal reports of any skin reactions, audited to ISO 9001:2015 by SGS standards.

- Hypoallergenic adhesive.
- Soft and flexible for a comfortable and secure wear.
- Stoma cap.
- Flat and convex high output pouches.
- Now available in drainable.

Flexible, hypoallergenic and skin friendly.





Flexible, hypoallergenic and skin friendly.



### Hydrocolloid Stoma Cap with Soft Cover and Filter

1420K Filter, opaque, soft cover both sides.

Pack Size: 50

Hole Size	Soft Cover - Opaque
Cut to fit 20 to 45mm	1420K

Flexible, hypoallergenic and skin friendly.



### Drainable Standard Hydrocolloid High Output Soft Convex

SCHV25-4115K - SCHV45-4115K Window, opaque, filter, soft cover both sides, wide bung closure.

Pack Size: 10

Hole Size	Plateau Size	Window
Cut to fit 15 to 24mm	25mm	SCHV25-4115K
Cut to fit 15 to 34mm	35mm	SCHV35-4115K
Cut to fit 15 to 44mm	45mm	SCHV45-4115K



### Drainable Standard Hydrocolloid High Output

HVP20-60 Window, opaque, filter, soft cover both sides, wide bung closure.

Pack Size: 10

Hole Size	Window
Cut to fit 20 to 60mm	HVP20-60



**WOUND  
MANAGEMENT**

Designed for Healthcare Professionals, our comprehensive range of Wound Management products includes Wound Managers and Option Managers.



Reliable wound management for Healthcare Professionals.



### Option Manager

OM-S Clear, soft cover one side, cap closure.  
Approx capacity 50ml.

OM-M Clear, soft cover one side, cap closure.  
Approx capacity 150ml.

OM-L Clear, soft cover one side, cap closure.  
Approx capacity 300ml.

Pack Size: 30

Hole Size	Window
Cuts up to 30mm	OM-S
Cuts up to 50mm	OM-M
Cuts up to 50mm	OM-L

Reliable wound management for Healthcare Professionals.



### Option Wound Manager Extra Small

WWM-XS - WWM-XSA Hole 10mm, film: clear.

Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 90mm wide x 180mm high	WWM-XS	WWM-XSA

(Oval shape 90mm x 180mm at maximum point)



### Option Wound Manager Small

WWM-S - WWM-SA Hole 10mm, film: clear.

Pack Size: 10

Hole Size	Standard	Access Port
Suitable for horizontal wounds up to 240mm wide x 155mm high	WWM-S	WWM-SA



### Option Wound Manager Oblong

WWM-O - WWM-OA Hole 10mm, film: clear.

Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 60mm wide x 190mm high	WWM-O	WWM-OA



### Option Wound Manager Medium

WWM-M - WWM-MA Hole 10mm, film: clear.

Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 85mm wide x 250mm high	WWM-M	WWM-MA

Reliable wound management for Healthcare Professionals.



### Option Wound Manager Square

WWM-SQ - WWM-SQA Hole 10mm, film: clear.

Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 160mm wide x 200mm high	WWM-SQ	WWM-SQA



### Option Wound Manager Large

WWM-L - WWM-LA Hole 10mm, film: clear.

Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 160mm wide x 260mm high	WWM-L	WWM-LA



### Option Wound Manager Extra Large\*

WWM-XL - WWM-XLA Hole 10mm, film: clear.

Pack Size: 5

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 230mm wide x 380mm high	WWM-XL	WWM-XLA

\*Option Wound Manager Extra Large is not available on prescription.



For added security and peace of mind the Oakmed comprehensive range of stoma care and related accessories includes solidifying capsules, tablets, gentle adhesive remover sprays and wipes, flange extenders and ostomy belts, and includes our new range of Trio Sil2™ accessories.

**NEW  
for  
2025**

Trio Sil2™ accessories are now available in our trusted accessories range.



Comprehensive range of stoma care accessories.

The Oakmed range of flange extenders now includes Genii™ silicone flange extenders using Sil2™ technology.



### Curvies

Secure and flexible silicone adhesive flange extenders. 110mm max width.

CUR001 Pack Size: 30

	Code
Curvies	CUR001



### Curvies Plus

Secure and flexible silicone adhesive, large flange extenders. 155mm max width.

CUR005 Pack Size: 20

	Code
Curvies Plus	CUR005



### Curvies Hydro

Kind to your skin, Hydrocolloid adhesive, large flange extenders. 155mm max width.

HYD007 Pack Size: 20

	Code
Curvies Hydro	HYD007



### Genii™ Silicone Flange Extenders

Provides a soft, flexible, invisible border which conforms to the body

TR1060 Pack Size: 20

	Code
Genii™ Flange Extenders	TR1060



Comprehensive range of stoma care accessories.



### EasyLease Sting Free Adhesive Remover Spray

Readily available in spray, EasyLease is the no sting medical adhesive remover with antiseptic tea tree oil. EasyLease gently helps to remove bags and clean any adhesive residue from the skin.

	Code
50ml Adhesive Remover Aerosol	1408



### Elite® Sting Free Adhesive Remover

100% silicone formulation. Hypoallergenic and sting free spray, non-irritant adhesive remover. Also available in a wipe.

	Code
50ml spray can	TR201
Pack of 30 wipes	TR202



### Elisse® Sting Free Skin Barrier

Hypoallergenic silicone formula that builds a protective barrier between the ostomy wafer or base plate and the skin. Also available in a wipe.

	Code
50ml spray can	TR204
Pack of 30 wipes	TR203



Comprehensive range of stoma care accessories.



### Option Seals

High tack and absorbent.

Pack of 10

	Code
Slim Seals - Slim and stretchy, ideal for skin protection.	WM-500
Seals - Slightly thicker than the slim seals for added depth.	WM-100



### Genii™ Ostomy Seals

Protects the skin immediately around the stoma from irritation and breakdown

Pack Size: 10

	Code
Small - 20mm-35mm	TR1020
Medium - 35mm-50mm	TR1028
Large - 50mm-80mm	TR1035



### Genii™ Convex Ostomy Seals

Soft and mildly convex to provide the best fit to the body shape. Can be applied convex side up or down.

Pack Size: 10

	Code
Small - 20mm-35mm	TR1020cx
Large - 30mm-50mm	TR1030cx





BROUGHT TO YOU BY

---

# OakMed

Giving you more choice.

[www.oakmed.co.uk](http://www.oakmed.co.uk)

**Oakmed Limited,**  
Unit B, The I O Centre, Barn Way,  
Lodge Farm Industrial Estate,  
Northampton, NN5 7UW  
Freephone: 0800 592 786  
Main: 01604 586 529 | Fax: 01604 584 784



Stay connected

Like us on facebook @OakmedUK

## GoldCare®

Home delivery by **OakMed**

Reliable stoma care delivery made easy

Freephone: 0800 592 786  
Register Online [goldcare.healthcare/register](http://goldcare.healthcare/register)

