

PRODUCT GUIDE



www.oakmed.co.uk



Oakmed Labels Explained

Please use this chart to help you identify individual products within the Oakmed range.

- All Oakmed products are colour coded and have a colour coded logo.
- You will see these colours and logos throughout our labelling, product guides, marketing and website.
- Some of the labels may have two colours if they come into two categories, e.g. Soft Convex^{*}/Alginate.
- All stock sent into hospitals are half stock boxes and state that on the label.

NB. Your Oakmed Representative is happy to check Oakmed stock for you when they visit.







Freephone: 0800 592 786 | oakmed.co.uk

New For 2022

Look out for

In our latest product development, we have updated our drainable range with our unique anatomical shape and a new smooth tab on our soft closure. Designed with patients and Healthcare Professionals to increase comfortable weartime.

Primo - Page 5 Soft Convex[®] - Page 8 Alginate - Page 18



New Alginate Wafers

Our healing Alginate wafers are now available in our Soft Convex and Low Profile Soft Convex, closed and drainable ranges.

Pages 18-20



New to our Seals range, Slim Seals are a cost effective way to manage a difficult stoma.



Freephone: 0800 592 786 | oakmed.co.uk

Contents

PRIMO®

Primo®

The perfect Hydrocolloid tapered edge for outstanding adhesion and security.

Page 5



Alginate

An adhesive which contains Alginate. Promotes skin integrity, absorbs moisture, and helps to heal sore skin.

Page 18



Hydrocolloid

A soft, flexible, hypoallergenic adhesive.

Page 33



Soft Convex®

7mm convexity. Available in tapered edge Hydrocolloid and Alginate adhesives. Also high output.

Page 8

Urostomy

Page 21

WOUND

Available in tap and soft bung

with Hydrocolloid, Microskin,

Soft Convex and Low Profile

Soft Convex adhesives.



Low Profile Soft Convex®

3.5mm convexity. Available in tapered edge Hydrocolloid and Alginate adhesives.

Page 12



Microskin®

The lightest, thinnest and most flexible wafer available.

Page 31



Wound Care

Wound and fistula drainage appliances in a range of sizes.

Page 36



Accessories

A comprehensive range of stoma and related accessories.

Page 40

PRIMO

Designed by patients and Healthcare Professionals, the unique features and benefits of our Primo[®] bags offer you the perfect fit with outstanding adhesion and assured security.

These features are integral to both our drainable and closed bags.

- Oakmed's Hydrocolloid has a long term reputation for minimal reports of any skin reactions, audited to ISO 9001:2015 by SGS.
- Anatomically shaped
- Infinitely superior, high tack adhesion.
- Leaves little or no residue when removed.
- Unique tapered edge, integral border.
- Hypoallergenic, non-irritant Hydrocolloid.
- Large adhesive area even when cut to maximum.
- Carbon filters designed individually for closed and drainable bag requirements.
- Improved soft touch closure on drainable bags.
- Soft non-cling fabric cover with viewing window.

Our most popular range with mini, midi, standard and maxi sizes available. Soft touch, easy to use closure on drainable bags.







Closed Mini Hydrocolloid Primo

HP11C10 Window, opaque, filter, soft cover both sides. Pack Size: 30

Hole Size	Window
Cut to fit 10 to 50mm	HP11C10



HP5C20 - HP5C60 Window, opaque, filter, soft cover both sides. HP5C20C Clear, filter, soft cover one side. Pack Size: 30

Hole Size	Window	Clear
Cut to fit 20 to 60mm	HP5C20	HP5C20C
Pre-cut to 25mm	HP5C25	
Pre-cut to 30mm	HP5C30	
Pre-cut to 35mm	HP5C35	
Pre-cut to 40mm	HP5C40	
Pre-cut to 45mm	HP5C45	
Pre-cut to 50mm	HP5C50	
Pre-cut to 55mm	HP5C55	
Pre-cut to 60mm	HP5C60	



Closed Maxi Hydrocolloid Primo

HP8C20 Window, opaque, filter, soft cover both sides. HP7C20C Clear, filter, soft cover one side. Pack Size: 30

Hole Size	Window	Clear
Cut to fit 20 to 90mm	HP8C20	HP7C20C

PRIMO'S U



Drainable Midi Hydrocolloid Primo

 $\mathsf{HP8D20}$ - $\mathsf{HP8D60}$ Window, opaque, filter, soft cover both sides, soft touch closure.

HP8D20C Clear, filter, soft cover one side, soft touch closure. Pack Size: 30

Hole Size	Window	Clear
Cut to fit 20 to 60mm	HP8D20	HP8D20C
Pre-cut to 25mm	HP8D25	
Pre-cut to 30mm	HP8D30	
Pre-cut to 35mm	HP8D35	
Pre-cut to 40mm	HP8D40	
Pre-cut to 45mm	HP8D45	
Pre-cut to 50mm	HP8D50	
Pre-cut to 55mm	HP8D55	
Pre-cut to 60mm	HP8D60	



Drainable Standard Hydrocolloid Primo

HP10D20 - HP10D60 Window, opaque, filter, soft cover both sides, soft touch closure.

HP10D20C Clear, filter, soft cover one side, soft touch closure. Pack Size: 30

Hole Size	Window	Clear
Cut to fit 10 to 50mm	HP10D20	HP10D20C
Pre-cut to 25mm	HP10D25	
Pre-cut to 30mm	HP10D30	
Pre-cut to 35mm	HP10D35	
Pre-cut to 40mm	HP10D40	
Pre-cut to 45mm	HP10D45	
Pre-cut to 50mm	HP10D50	
Pre-cut to 55mm	HP10D55	
Pre-cut to 60mm	HP10D60	

Soft Convex®

7mm convexity

Naturally the flexible choice for maximum security and added confidence, our Soft Convex[®] range of bags fit the contours of your body for more comfortable wear. In both drainable and closed, tapered edge Hydrocolloid and healing Alginate adhesives, our range also offers a high output bag.

Unique features:

- 7mm depth of convexity.
- Maintains its convexity whilst
 in use.
- A comprehensive range of cut to fit and pre-cut sizes.
- High tack adhesion for advanced performance.
- Skin friendly design, hypoallergenic adhesives.
- Efficient odour eliminating filter system.

NEW FOR 2022


Unique anatomical shape now on both closed and drainable bags. Soft touch, easy to use closure on drainable bags.





Naturally the flexible choice.





Closed Standard Hydrocolloid

Soft Convex



Drainable Midi Hydrocolloid Soft Convex

 $\mathsf{SC25-4315KV}$ - $\mathsf{SC45-4340KV}$ Window, opaque, filter, soft cover both sides, soft touch closure.

 $\mathsf{SC25}\text{-}4315\mathsf{KVC}$ - $\mathsf{SC45}\text{-}4315\mathsf{KVC}$ Clear, filter, soft cover one side, soft touch closure.

Pack Size: 20

Hole Size	Plateau Size	Window
Cut to fit 15 to 24mm	25mm	SC25-4315KV
Pre-cut to 20mm	25mm	SC25-4320KV
Cut to fit 15 to 34mm	35mm	SC35-4315KV
Pre-cut to 25mm	35mm	SC35-4325KV
Pre-cut to 30mm	35mm	SC35-4330KV
Cut to fit 15 to 44mm	45mm	SC45-4315KV
Pre-cut to 35mm	45mm	SC45-4335KV
Pre-cut to 40mm	45mm	SC45-4340KV
Hole Size	Plateau Size	Clear

	Plateau Size	Clear
Cut to fit 15 to 24mm	25mm	SC25-4315KVC
Cut to fit 15 to 34mm	35mm	SC35-4315KVC
Cut to fit 15 to 44mm	45mm	SC45-4315KVC

7mm depth





Drainable Standard Hydrocolloid Soft Convex

 $\mathsf{SC25-4120KV}$ - $\mathsf{SC57-4120KV}$ Window, opaque, filter, soft cover both sides, soft touch closure.

Pack Size: 20

Hole Size	Plateau Size	Window
Cut to fit 20 to 24mm	25mm	SC25-4120KV
Cut to fit 20 to 34mm	35mm	SC35-4120KV
Pre-cut to 25mm	35mm	SC35-4125KV
Pre-cut to 30mm	35mm	SC35-4130KV
Cut to fit 20 to 44mm	45mm	SC45-4120KV
Pre-cut to 35mm	45mm	SC45-4135KV
Pre-cut to 40mm	45mm	SC45-4140KV
Pack Size: 10		
Cut to fit 20 to 57mm	57mm	SC57-4120KV



Naturally the flexible choice.



Closed Standard Alginate Soft Convex

 $\mathsf{SCA25}\text{-}\mathsf{O420K}$ - $\mathsf{SCA45}\text{-}\mathsf{O440K}$ Window, opaque filter, soft cover both sides.

Pack Size: 20

Plateau Size	Window
25mm	SCA25-0420K
35mm	SCA35-0420K
35mm	SCA35-0425K
35mm	SCA35-0430K
45mm	SCA45-0420K
45mm	SCA45-0435K
45mm	SCA45-0440K
	25mm 35mm 35mm 35mm 45mm



Drainable Standard Hydrocolloid High Output Soft Convex

SCHV25-4115K - SCHV45-4115K Window, opaque, filter, soft cover both sides, wide bung closure.

Pack Size: 10

Hole Size	Plateau Size	Window
Cut to fit 15 to 24mm	25mm	SCHV25-4115K
Cut to fit 15 to 34mm	35mm	SCHV35-4115K
Cut to fit 15 to 44mm	45mm	SCHV45-4115K



Drainable Standard Alginate Soft Convex

 $\mathsf{SCA25}\xspace{-}4120\mathsf{KV}$ - $\mathsf{SCA45}\xspace{-}4140\mathsf{KV}$ Window, opaque, filter, soft cover both sides, soft touch closure.

Pack Size: 20

Hole Size	Plateau Size	Window
Cut to fit 20 to 24mm	25mm	SCA25-4120KV
Cut to fit 20 to 34mm	35mm	SCA35-4120KV
Pre-cut to 25mm	35mm	SCA35-4125KV
Pre-cut to 30mm	35mm	SCA35-4130KV
Cut to fit 20 to 44mm	45mm	SCA45-4120KV
Pre-cut to 35mm	45mm	SCA45-4135KV
Pre-cut to 40mm	45mm	SCA45-4140KV

7mm depth



Low Profile Soft Convex

3.5mm convexity

Offering the same superior flexibility and comfortable wear, our Low Profile Soft Convex[®] range gives you 3.5mm of convexity.

Available in both drainable and closed, tapered edge Hydrocolloid and healing Alginate adhesives.

Unique features:

- 3.5mm depth of convexity.
- Extremely cost effective product.
- Soft and flexible design, works like a flat bag with a seal.
- A comprehensive range of cut to fit and pre-cut sizes.
- High tack adhesion for advanced performance.
- Skin friendly design, hypoallergenic adhesives.
- Efficient odour eliminating filter system.

NEW FOR 2022

Unique anatomical shape now on closed and drainable bags, designed in conjunction with Healthcare Professionals.







Closed Standard Hydrocolloid Low Profile Soft Convex

LPSC28-0420K - LPSC47-0445K Window, opaque, filter, soft cover both sides. Pack Size: 20

Plateau Size	Window
28mm	LPSC28-0420K
28mm	LPSC28-0425K
40mm	LPSC40-0420K
40mm	LPSC40-0430K
40mm	LPSC40-0435K
47mm	LPSC47-0420K
47mm	LPSC47-0440K
47mm	LPSC47-0445K
	28mm 28mm 40mm 40mm 40mm 47mm

3.5mm depth



Low Profile Soft Convex



Drainable Midi Hydrocolloid Low Profile Soft Convex

LPSC28-4315KV - LPSC47-4345KV Window, opaque, filter, soft cover both sides, soft touch closure. LPSC28-4315KVC - LPSC47-4315KVC Clear, filter, soft cover one side, soft touch closure. Pack Size: 20

Hole Size	Plateau Size	Window	Clear
Cut to fit 15 to 28mm	28mm	LPSC28-4315KV	LPSC28-4315KVC
Pre-cut to 20mm	28mm	LPSC28-4320KV	
Pre-cut to 25mm	28mm	LPSC28-4325KV	
Cut to fit 15 to 40mm	40mm	LPSC40-4315KV	LPSC40-4315KVC
Pre-cut to 30mm	40mm	LPSC40-4330KV	
Pre-cut to 35mm	40mm	LPSC40-4335KV	
Cut to fit 15 to 47mm	47mm	LPSC47-4315KV	LPSC47-4315KVC
Pre-cut to 40mm	47mm	LPSC47-4340KV	
Pre-cut to 45mm	47mm	LPSC47-4345KV	







Drainable Standard Hydrocolloid Low Profile Soft Convex

LPSC28-4115KV - LPSC47-4145KV Window, opaque, filter, soft cover both sides, soft touch closure. LPSC28-4115KVC - LPSC47-4115KVC Clear, filter, soft cover one side, soft touch closure. Pack Size: 20

Hole Size	Plateau Size	Window	Clear
Cut to fit 15 to 28mm	28mm	LPSC28-4115KV	LPSC28-4115KVC
Pre-cut to 20mm	28mm	LPSC28-4120KV	
Pre-cut to 25mm	28mm	LPSC28-4125KV	
Cut to fit 15 to 40mm	40mm	LPSC40-4115KV	LPSC40-4115KVC
Pre-cut to 30mm	40mm	LPSC40-4130KV	
Pre-cut to 35mm	40mm	LPSC40-4135KV	
Cut to fit 15 to 47mm	47mm	LPSC47-4115KV	LPSC47-4115KVC
Pre-cut to 40mm	47mm	LPSC47-4140KV	
Pre-cut to 45mm	47mm	LPSC47-4145KV	



Low Profile Soft Convex



3.5mm depth

Closed Standard Alginate Low Profile Soft Convex

LPSCA28-0420K - LPSCA47-0445K Window, opaque, filter, soft cover both sides. Pack Size: 20

Hole Size	Plateau Size	Window
Cut to fit 20 to 28mm	28mm	LPSCA28-0420K
Pre-cut to 25mm	28mm	LPSCA28-0425K
Cut to fit 20 to 40mm	40mm	LPSCA40-0420K
Pre-cut to 30mm	40mm	LPSCA40-0430K
Pre-cut to 35mm	40mm	LPSCA40-0435K
Cut to fit 20 to 47mm	47mm	LPSCA47-0420K
Pre-cut to 40mm	47mm	LPSCA47-0440K
Pre-cut to 45mm	47mm	LPSCA47-0445K



Naturally the flexible choice.





Drainable Standard Alginate Low Profile Soft Convex

LPSCA28-4115KV - LPSCA47-4145KV Window, opaque, filter, soft cover both sides, soft touch closure. Pack Size: 20

Hole Size	Plateau Size	Window
Cut to fit 15 to 28mm	28mm	LPSCA28-4115KV
Pre-cut to 20mm	28mm	LPSCA28-4120KV
Pre-cut to 25mm	28mm	LPSCA28-4125KV
Cut to fit 15 to 40mm	40mm	LPSCA40-4115KV
Pre-cut to 30mm	40mm	LPSCA40-4130KV
Pre-cut to 35mm	40mm	LPSCA40-4135KV
Cut to fit 15 to 47mm	47mm	LPSCA47-4115KV
Pre-cut to 40mm	47mm	LPSCA47-4140KV
Pre-cut to 45mm	47mm	LPSCA47-4145KV



Our Alginate wafer is one of our superior skin protectors, offering both preventative and curative properties. Available in a wide range of sizes in both drainable and closed bags.

Used in woundcare for many years. Alginates are:

• Healing.

Alginate

- Absorbent.
- Anti-bacterial.
- Haemostatic.
- Used without any accessories for a cost effective way to manage problematic skin.

Application for use are:

- Skin prone to allergies.
- Sore, excoriated peri-stomal skin.
- Moist bleeding, peri-stomal skin.
- Ulcerated peri-stomal skin.
- Muco-cutaneous separation.
- Stomas which produce excessive mucous.
- Infected peri-stomal suture sites.
- Alginate adhesive can readily be used on non-problematic skin to prevent complications occurring.

Oakmed Alginate adhesives contain not just a trace element but **15%** Alginate, this is significant enough to promote healing.





With 15% Alginate, stimulates healing for sore skin.



Hole Size	Window
Cut to fit 10 to 50mm	AHP11C10

AHP11C10 Window, opaque, filter, soft cover both sides.



Closed Standard Alginate

Closed Mini Alginate

Pack Size: 30

AHP5C20 - AHP5C60 Window, opaque, filter, soft cover both sides. Pack Size: 30

Hole Size	Window
Cut to fit 20 to 60mm	AHP5C20
Pre-cut to 25mm	AHP5C25
Pre-cut to 30mm	AHP5C30
Pre-cut to 35mm	AHP5C35
Pre-cut to 40mm	AHP5C40
Pre-cut to 45mm	AHP5C45
Pre-cut to 50mm	AHP5C50
Pre-cut to 55mm	AHP5C55
Pre-cut to 60mm	AHP5C60



Closed Maxi Alginate

AHP8C20 Window, opaque, filter, soft cover both sides. AHP7C20C Clear, soft cover one side. Pack Size: 30

Hole Size	Window	Clear
Cut to fit 20 to 90mm	AHP8C20	AHP7C20C

With **15% Alginate**, stimulates healing for sore skin.



Alginate

Drainable Midi Alginate

AV-4320K Window, filter, soft cover both sides, soft touch closure.

Pack Size: 30

Hole Size	Window Midi
Cut to fit 20 to 60mm	AV-4320K



Drainable Standard Alginate

AV-4120K - AV-4160K Window, filter, soft cover both sides, soft touch closure. Pack Size: 30

Hole Size	Window
Cut to fit 20 to 60mm	AV-4120K
Pre-cut to 25mm	AV-4125K
Pre-cut to 30mm	AV-4130K
Pre-cut to 35mm	AV-4135K
Pre-cut to 40mm	AV-4140K
Pre-cut to 45mm	AV-4145K
Pre-cut to 50mm	AV-4150K
Pre-cut to 55mm	AV-4155K
Pre-cut to 60mm	AV-4160K



Our extended Urostomy range is shaped for maximum capacity with minimum bulk. The unique anti-twist shape allows for easier emptying. Now available in soft bung or new tap closure for added confidence and convenience.

- Available in a choice of adhesives.
- Soft fabric cover, quiet, non-cling, comfortable.
- Split fabric cover for easy positioning.
- Non-return valve, stops urine back flowing.
- Soft bung closure, safe from leaks.
- Anti-twist shape.
- Drip cap on tap closure.

NEW TAP CLOSURE FOR 2022

Ask about our Clinic/Community packs, the easy way for patients to find the right bag for them.





Maximum capacity, minimum bulk for increased confidence and comfort.





Maximum capacity, minimum bulk for increased confidence and comfort.





Urostomy Standard Hydrocolloid

GC1500 Window, opaque, soft cover both sides, soft bung closure. Pack Size: 20

	Hole Size	Window
*	Cut to fit 15 to 60mm	GC1500



Hole Size			Window
Pack Size:	20		
Pack Size:	20		
9013001	willow, opaque, sont cover	UUTIT SIDES	, tap closure

	Hole Size	window
k	Cut to fit 15 to 60mm	GC1500T

Urostomy Standard Hydrocolloid

GC1500C Clear, soft cover one side, soft bung closure. Pack Size: 20

	Hole Size	Clear
*	Cut to fit 15 to 60mm	GC1500C



Available as a Clinic/Community pack. Call us to order.



Maximum capacity, minimum bulk, with 7mm convexity, for increased confidence and comfort.



SCU25-0015

SCU25-0020

SCU35-0015

SCU35-0025

SCU35-0030

SCU45-0015

SCU45-0035

SCU45-0040



Maximum capacity, minimum bulk, with 7mm convexity, for increased confidence and comfort.





Urostomy Standard Hydrocolloid Soft Convex

SCU25-0015T - SCU45-0040T Window, opaque, soft cover

		Hole Size	Plateau Size	Window
11	*	Cut to fit 15 to 24mm	25mm	SCU25-0015T
		Pre-cut to 20mm	25mm	SCU25-0020T
	*	Cut to fit 15 to 34mm	35mm	SCU35-0015T
1		Pre-cut to 25mm	35mm	SCU35-0025T
1		Pre-cut to 30mm	35mm	SCU35-0030T
	*	Cut to fit 15 to 44mm	45mm	SCU45-0015T
		Pre-cut to 35mm	45mm	SCU45-0035T
		Pre-cut to 40mm	45mm	SCU45-0040T



Urostomy

Urostomy Standard Hydrocolloid Soft Convex

Urostomy Standard Hydrocolloid Soft Convex

Plateau Size

25mm

25mm

35mm

35mm

35mm

45mm

45mm

45mm

SCU25-0015 - SCU45-0040 Window, opaque, soft cover

both sides, soft bung closure.

Pack Size: 20

Hole Size 📯 Cut to fit 15 to 24mm

Pre-cut to 20mm

📯 Cut to fit 15 to 34mm

Pre-cut to 25mm

Pre-cut to 30mm

Pre-cut to 35mm

Pre-cut to 40mm

Cut to fit 15 to 44mm

R

SCU25-0015C - SCU45-0015C Clear, soft cover one side, soft bung closure. Pack Size: 20

	Hole Size	Plateau Size	Clear
*	Cut to fit 15 to 24mm	25mm	SCU25-0015C
*	Cut to fit 15 to 34mm	35mm	SCU35-0015C
*	Cut to fit 15 to 44mm	45mm	SCU45-0015C



Urostomy Standard Hydrocolloid Soft Convex

SCU25-0015CT - SCU45-0015CT Clear, soft cover one side, tap closure.

Pack Size: 20

	Hole Size	Plateau Size	Clear
*	Cut to fit 15 to 24mm	25mm	SCU25-0015CT
*	Cut to fit 15 to 34mm	35mm	SCU35-0015CT
*	Cut to fit 15 to 44mm	45mm	SCU45-0015CT

Available as a Clinic/Community pack. Call us to order.



Maximum capacity, minimum bulk, with 3.5mm convexity, for increased confidence and comfort.



Urostomy Standard Hydrocolloid Low Profile Soft Convex

LPSCU28-0015 - LPSCU47-0045 Window, opaque, soft cover both sides, bung closure.

Pack Size: 20

	Hole Size	Plateau Size	Window
*	Cut to fit 15 to 28mm	28mm	LPSCU28-0015
	Pre-cut to 20mm	28mm	LPSCU28-0020
	Pre-cut to 25mm	28mm	LPSCU28-0025
*	Cut to fit 15 to 40mm	40mm	LPSCU40-0015
	Pre-cut to 30mm	40mm	LPSCU40-0030
	Pre-cut to 35mm	40mm	LPSCU40-0035
*	Cut to fit 15 to 47mm	47mm	LPSCU47-0015
	Pre-cut to 40mm	47mm	LPSCU47-0040
	Pre-cut to 45mm	47mm	LPSCU47-0045



Urostomy Standard Hydrocolloid Low Profile Soft Convex

LPSCU28-0015C - LPSCU47-0015C Clear, soft cover one side, soft bung closure. Pack Size: 20

	Hole Size	Plateau Size	Clear
*	Cut to fit 15 to 28mm	28mm	LPSCU28-0015C
*	Cut to fit 15 to 40mm	40mm	LPSCU40-0015C
*	Cut to fit 15 to 47mm	47mm	LPSCU47-0015C

Available as a Clinic/Community pack. Call us to order.



Maximum capacity, minimum bulk, with 3.5mm convexity, for increased confidence and comfort.

NEW



Urostomy Standard Hydrocolloid Low Profile Soft Convex

LPSCU28-0015T - LPSCU47-0045T Window, opaque, soft cover both sides, tap closure.

Pack Size: 20

	Hole Size	Plateau Size	Window
*	Cut to fit 15 to 28mm	28mm	LPSCU28-0015T
	Pre-cut to 20mm	28mm	LPSCU28-0020T
	Pre-cut to 25mm	28mm	LPSCU28-0025T
*	Cut to fit 15 to 40mm	40mm	LPSCU40-0015T
	Pre-cut to 30mm	40mm	LPSCU40-0030T
	Pre-cut to 35mm	40mm	LPSCU40-0035T
*	Cut to fit 15 to 47mm	47mm	LPSCU47-0015T
	Pre-cut to 40mm	47mm	LPSCU47-0040T
	Pre-cut to 45mm	47mm	LPSCU47-0045T



Urostomy Standard Hydrocolloid Low Profile

LPSCU28-0015CT - LPSCU47-0015CT Clear. soft cover

	Hole Size	Plateau Size	Clear
*	Cut to fit 15 to 28mm	28mm	LPSCU28-0015CT
*	Cut to fit 15 to 40mm	40mm	LPSCU40-0015CT
*	Cut to fit 15 to 47mm	47mm	LPSCU47-0015CT



Starter Pack

Clinic/community

Clinic/Community starter packs in all cut to fit Urostomy products

We have a full range of clinic starter packs available in both new 'Twist and Safe' tap and soft bung closure options.

Ideal to help Urostomists choose the right bag for their needs. Healthcare Professionals can rest assured that the pack contains everything the patient needs to trial our bag options. Full instructions for use are included to provide guidance on how to use the connector with the night drainage bag.

Please see over the page for the options we offer.

- 5 x Bags of chosen product
- I x Night drainage bag
- 2 x Connectors
- I x Belt
- I × IFU
- | x Product code information sheet


Ask your Oakmed Territory Manager for more details.

Please contact our customer service department on 0800 592 786 if you have any questions or need any help at all regarding this or any other Oakmed products.

		Hole Size	
	*	Cut to fit 15 to 60mm	
	*	Cut to fit 15 to 60mm	
TAP	*	Cut to fit 15 to 24mm	SC
	*	Cut to fit 15 to 34mm	SC
	*	Cut to fit 15 to 44mm	SC
	*	Cut to fit 15 to 24mm	
DRIP CAP	*	Cut to fit 15 to 34mm	
CAF	*	Cut to fit 15 to 44mm	
	*	Cut to fit 15 to 28mm	LPS
	*	Cut to fit 15 to 40mm	LPS
	*	Cut to fit 15 to 47mm	LPS
	*	Cut to fit 15 to 28mm	
	*	Cut to fit 15 to 40mm	
	*	Cut to fit 15 to 47mm	

Bung

Tap





	Hole Size	Window	Clear
	Cut to fit 15 to 60mm	GC1500	-
	Cut to fit 15 to 60mm	-	GC1500C
*	Cut to fit 15 to 24mm	SCU25-0015	-
*	Cut to fit 15 to 34mm	SCU35-0015	-
*	Cut to fit 15 to 44mm	SCU45-0015	-
*	Cut to fit 15 to 24mm	-	SCU25-0015C
*	Cut to fit 15 to 34mm	-	SCU35-0015C
*	Cut to fit 15 to 44mm	-	SCU45-0015C
*	Cut to fit 15 to 28mm	LPSCU28-0015	-
*	Cut to fit 15 to 40mm	LPSCU40-0015	-
*	Cut to fit 15 to 47mm	LPSCU47-0015	-
*	Cut to fit 15 to 28mm	-	LPSCU28-0015C
*	Cut to fit 15 to 40mm	-	LPSCU40-0015C
*	Cut to fit 15 to 47mm	-	LPSCU47-0015C



Maximum capacity, minimum bulk for increased confidence and comfort.





Urostomy Standard Microskin

MJ1300T - MJ3800T Window, opaque, soft cover both sides, tap closure.

Pack Size: 20

Hole Size	Window
Cut to fit 10 to 50mm	MJ1300T
Pre-cut to 16mm	MJ1600T
Pre-cut to 19mm	MJ1900T
Pre-cut to 22mm	MJ2200T
Pre-cut to 25mm	MJ2500T
Pre-cut to 32mm	MJ3200T
Pre-cut to 38mm	MJ3800T



Urostomy Standard Microskin

MJ1300 - MJ3800 Window, opaque, soft cover both sides, soft bung closure.

Pack Size: 20

Hole Size	Window
Cut to fit 10 to 50mm	MJ1300
Pre-cut to 16mm	MJ1600
Pre-cut to 19mm	MJ1900
Pre-cut to 22mm	MJ2200
Pre-cut to 25mm	MJ2500
Pre-cut to 32mm	MJ3200
Pre-cut to 38mm	MJ3800

Micro[®] Skin

Unique to Oakmed, the Microskin[®] adhesive is thin, transparent and flexible, moulding easily to all body contours and creases, significantly reducing the risk of leakage.

Unique features:

- Transparent, thin and flexible design.
- Moulds easily and comfortably to body contours.
- Barely visible on skin when in place.
- Breathable and moisture vapour permeable for healthier skin.
- Durable for extended wear time.

Transparent and flexible, Microskin® is the ideal wafer for active people.





Transparent and flexible, reduces the risk of leakage.



Closed Standard Microskin

M-0320K - M-0340K Window, filter, soft cover both sides. Pack Size: 30

Window
M-0320K
M-0325K
M-0330K
M-0335K
M-0340K



Closed Plus Microskin

M-0420K - M-0440K Window, filter, soft cover both sides. Pack Size: 30

Hole Size	Window
Cut to fit 20 to 50mm	M-0420K
Pre-cut to 25mm	M-0425K
Pre-cut to 30mm	M-0430K
Pre-cut to 35mm	M-0435K
Pre-cut to 40mm	M-0440K



Oakmed's Hydrocolloid has a long term reputation for minimal reports of any skin reactions, audited to ISO 9001:2015 by SGS standards.

- Hypoallergenic adhesive.
- Soft and flexible for a comfortable and secure wear.
- Stoma cap.
- Flat and convex high output pouches.
- Now available in drainable.



Drainable Standard Microskin

MV-4120K Window, filter, soft cover both sides, soft touch closure. Pack Size: 30

Hole Size	Standard
Cut to fit 20 to 50mm	MV-4120K







32



Flexible, hypoallergenic and skin friendly.



Hydrocolloid Stoma Cap with Soft Cover and Filter

1420K Filter, opaque, soft cover both sides. Pack Size: 50

Hole SizeSoft Cover - OpaqueCut to fit 20 to 50mm1420K



Drainable Standard Hydrocolloid

V-4120K - V-4160K Window, opaque, filter, soft cover both sides, soft touch closure.

Pack Size: 30

Window
V-4120K
V-4125K
V-4130K
V-4135K
V-4140K
V-4145K
V-4150K
V-4155K
V-4160K



Flexible, hypoallergenic and skin friendly.



Drainable Standard Hydrocolloid High Output

 $\operatorname{HVP20-60}$ Window, opaque, filter, soft cover both sides, wide bung closure.

Pack Size: 10

Hole Size	Window
Cut to fit 20 to 60mm	HVP20-60



Drainable Standard Hydrocolloid High Output Soft Convex

SCHV25-4115K - SCHV45-4115K Window, opaque, filter, soft cover both sides, wide bung closure.

Pack Size: 10

Hole Size	Plateau Size	Window
Cut to fit 15 to 24mm	25mm	SCHV25-4115K
Cut to fit 15 to 34mm	35mm	SCHV35-4115K
Cut to fit 15 to 44mm	45mm	SCHV45-4115K

34



Designed for Healthcare Professionals, our comprehensive range of Wound Management products includes Wound Managers and Option Managers.







Reliable wound management for Healthcare Professionals.



Option Manager

OM-S Clear, soft cover one side, cap closure. Approx capacity 50ml.

OM-M Clear, soft cover one side, cap closure. Approx capacity 150ml.

OM-L Clear, soft cover one side, cap closure. Approx capacity 300ml.

Pack Size: 30

Hole Size	Window
Cuts up to 30mm	OM-S
Cuts up to 50mm	OM-M
Cuts up to 50mm	OM-L

Reliable wound management for Healthcare Professionals.



Option Wound Manager Extra Small

WWM-XS - WWM-XSA Hole 10mm, film: clear. Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 90mm wide x 180mm high	WWM-XS	WWM-XSA

(Oval shape 90mm x 180mm at maximum point)



Option Wound Manager Small

WWM-S - WWM-SA Hole 10mm, film: clear. Pack Size: 10

Option Wound Manager Oblong WWM-O - WWM-OA Hole 10mm, film: clear.

Hole Size	Standard	Access Port
Suitable for horizontal wounds up to 240mm wide x 155mm high	WWM-S	WWM-SA



38

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 60mm wide x 190mm high	WWM-0	WWM-OA



WWM-M - WWM-MA Hole 10mm, film: clear. Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 85mm wide x 250mm high	WWM-M	WWM-MA



Reliable wound management for Healthcare Professionals.



Option Wound Manager Square

WWM-SQ - WWM-SQA Hole 10mm, film: clear. Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 160mm wide x 200mm high	WWM-SQ	WWM-SQA



Option Wound Manager Large

WWM-L - WWM-LA Hole 10mm, film: clear. Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 160mm wide x 260mm high	WWM-L	WWM-LA



Option Wound Manager Extra Large*

WWM-XL - WWM-XLA Hole 10mm, film: clear.

Pack Size: 5

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 230mm wide x 380mm high	WWM-XL	WWM-XLA

[•]Option Wound Manager Extra Large is not available on prescription.

Pack Size: 10



For added security and peace of mind the Oakmed comprehensive range of stoma care and related accessories includes solidifying capsules, tablets, gentle adhesive remover sprays and wipes, flange extenders and ostomy belts.

New Slim Seals and Hydrocolloid Curvies, larger and skin friendly for increased wear time.

40









Comprehensive range of stoma care accessories.

CURVIES

The Oakmed range of flange extenders now includes larger sizes, in our secure silicone and caring hypoallergenic Hydrocolloid adhesives.



Curvies

Secure and flexible silicone adhesive flange extenders. 110mm max width.

CUR001 Pack Size: 30

	Code
Curvies	CUR001





Secure and flexible silicone adhesive, large flange extenders. 155mm max width.

CUR005 Pack Size: 20

	Code
Curvies Plus	CUR005



Curvies Hydro

Kind to your skin, Hydrocolloid adhesive, large flange extenders. 155mm max width.

HYD007 Pack Size: 20

	Code
Curvies Hydro	HYD007

For latest product information visit our website **oakmed.co.uk**

Comprehensive range of stoma care accessories.



EasyLease[™]

Readily available in both spray and handy individual wipes, EasyLease is the no sting medical adhesive remover with antiseptic tea tree oil. EasyLease gently helps to remove bags and clean any adhesive residue from the skin.

	Code
Adhesive Remover Wipes	1206
Adhesive Remover Aerosol	1408



GelX

For use in Ostomy Bags. 140 average contents.

	Code
Solidifying Capsules	GEL-X
Solidifying Tablets	GEL-XT



Option Seals

Pack of 10

		Code
	Seals	WM-100
	Slim Seals	WM-500
_	Still Sedis	WIN-2000

Comprehensive range of stoma care accessories.



Ostomy Belt	
	Code
1 per pack	OCB02



Urostomy Drainage Adaptor

5

	Code
per pack	UDA



Available on prescription, ILex Skin Protectant is a topical skin barrier designed to protect the skin from harmful irritants such as urine and faeces.

- Proven 0.00cm/hour penetration rate. Forms a complete seal over sore skin and repels harmful fluids and bacteria.
- Soothes and heals damaged skin.
- Ideal for fistulas, ostomy site irritations and excoriation.

Website: www.ilexhealthproducts.com Telephone: +44 (0) 1604 645937



	Code
7g Tubes - Box of 20	IPH
57g Tube	IP5I*
227g Tub	IPT50

[•]Code IP51 (57g Tube) available on prescription.

BROUGHT TO YOU BY



www.oakmed.co.uk

Oakmed Limited, Unit B, The I O Centre, Barn Way, Lodge Farm Industrial Estate, Northampton, NN5 7UW

Freephone: 0800 592 786 Main: 01604 586 529 | Fax: 01604 584 784





Your fast, FREE & simple stoma care products delivery service

Freephone: 0800 592 786 Register Online goldcare.healthcare/register

